Waist and abdomen massager

Instructions

 Danger	
Do not use it at the same time as other medical electronic devices, such as pacemakers artificial hearts and lungs and other life-sustaining medical electronic devices, electrocardiographs and other medical electronic devices, otherwise it will cause danger. If the patient uses both high-frequency surgical equipment and massagers, burns may be caused at the massager blade position and may damage the massager: if the massager is used near the shortwave or microwave massage equipment (1 meter), the output of the massager may be unstable, and the use of electrodes close to the chest increases the risk of cardiac fibrillation.	
 Danger	
Pregnant women and women who have menstruation, sensitive skin, heart disease, abnormal blood pressure, malignant tumors, cerebrovascular patients, acute patients or other patients who are receiving medical treatment must consult a doctor before using this product, and do not use it if you are sensitive to electricity.	
1. Do not use near the heart, head, mouth or genitals, skin diseases and other parts. 2. Contraindicated in patients with impaired skin perception or in patients who are not sensitive to heat 3. This product cannot be used in high heat, flammable, electromagnetic radiation, humid places. 4. It is forbidden to use when bathing, sweating, and sleeping.	
Do not disassemble, repair and modify the massager without authorization it will cause malfunction or electric shock.	
 Be careful	
1. Do not move the massager during use, when changing the massage part, it should be turned off before use, otherwise there will be strong stimulation. 2. Do not use it on children or people without the ability to express consciousness. 3. If you feel unwell due to the use of this product, please stop using it immediately and consult a doctor. 4. Please turn off the power when not using the product. 5. Do not use this product beyond the scope.	

Operation and use method

002

004

1.002 How to turn on

- ✧ Plug in the Android mobile phone data cable. After power on, the instrument enters the standby state, and the red light indicator is always on.
- ✧ Long press the "switch" key on the instrument for 1 second to start the machine, the instrument starts to work, the blue indicator is always on, and turn on frequency conversion vibration massage and infrared physiotherapy.
- ✧ The intensity of frequency conversion vibration massage and infrared physiotherapy is divided into three grades, which will be automatically adjusted to the highest grade after startup. After startup, briefly press the "switch" key to adjust the gear.
- ✧ Press and hold the "switch" key on the instrument for 1 second to turn it off.

2.002 How to connect the power supply

- ✧ Please use the USB charging cable provided with the instrument to connect the mobile phone charger or charging treasure (the output voltage is 5V and the current is not greater than 1.5A). After power on, the red light is on and the blue light is on after work.

3.002 Why doesn't the product have a charging adapter

- ✧ The product supports universal Android USB data cable
- ✧ Support the standard universal mobile phone adapter and powerbank on the market.
(Note: Output Voltage 5V, output current > 500mA, recommended output current > = 1a)

1.004 How to turn on

- ✧ Plug in the Android mobile phone data cable. After power on, the instrument enters the standby state, and the red light indicator is always on.
- ✧ Press the "switch" key on the instrument to start up, the instrument starts to work, the blue indicator light is always on, starts vibration massage, left and right massage, frequency conversion vibration.
- ✧ Press the "switch" key on the instrument to shut down.
- ✧ During use, press the "+" key once to increase the strength by one gear, Each time you press the "-" key, the strength will be reduced by one gear. There are 3 gears in total. Please choose your own comfort intensity.

2.004 How to connect the power supply

- ✧ Please use the USB charging cable provided with the instrument to connect the mobile phone charger or charging treasure (the output voltage is 5V and the current is not greater than 1.0A). After power on, the red light of the indicator is on and the blue light is on after work.

3.004 Why doesn't the product have a charging adapter

- ✧ The product supports universal Android USB data cable
- ✧ Support the standard universal mobile phone adapter and powerbank on the market.
(Note: Output Voltage 5V, output current > 500mA, recommended output current > = 1a)

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